

# PLWS Specifications

Low Calf Support position:	16.5" from floor
Low Stirrup position:	20.5" from floor
High Calf Support position:	29" from floor
High Stirrup position:	33" from floor
Stirrup inside width:	5.75"
Calf support surface width:	8.75"
<hr/>	
Shipping Weight Product:	9.0 lbs (4.1 kg)
Shipping Weight Boxed:	10.4 lbs (4.7 kg)
<hr/>	
Tested to support 150 lbs.	

Listed to the Safety Standards for Medical Electrical Equipment UL 60601-1, CSA-C22.2 60601-1, and IEC 60601-1.

